


Indications

- Peripheral arterial reconstructions
- Extra-anatomical procedures
- Dialysis shunt placement

VascuGraft® NEO

- Unique VascuGraft® wall design with well-defined and different fibril lengths for a low thrombogenicity.
- Additional ePTFE wrap supports enhanced suture retention and high burst strength.
- Clear guideline indicates diameter and wall thickness of the graft.
- Helical reinforcement for high resistance against kinking and compression.
- Excellent and pliable handling and suture behaviour. Simple and easy removal of the spiral reinforcement.



	Length	Size	Article code	pcs.	
Standard wall	20 cm	5 mm	V1103505	1	
		6 mm	V1103506	1	
		7 mm	V1103507	1	
		8 mm	V1103508	1	
	30 cm	5 mm	V1103510	1	
	40 cm	6 mm	V1103516	1	
		7 mm	V1103517	1	
		8 mm	V1103518	1	
	50 cm	5 mm	V1103520	1	
	60 cm	6 mm	V1103526	1	
		7 mm	V1103527	1	
		8 mm	V1103528	1	
	70 cm	5 mm	V1103530	1	
	80 cm	6 mm	V1103536	1	
		7 mm	V1103537	1	
		8 mm	V1103538	1	
Standard wall, helix	30 cm	5 mm	V1103560	1	
	40 cm	6 mm	V1103566	1	
		7 mm	V1103567	1	
		8 mm	V1103568	1	
	50 cm	5 mm	V1103570	1	
	60 cm	6 mm	V1103576	1	
		7 mm	V1103577	1	
		8 mm	V1103578	1	
	70 cm	5 mm	V1103580	1	
	80 cm	6 mm	V1103586	1	
		7 mm	V1103587	1	
		8 mm	V1103588	1	
		10 mm	V1103589	1	
	Thin wall	30 cm	5 mm	V1103610	1
40 cm		6 mm	V1103616	1	
		7 mm	V1103617	1	
		8 mm	V1103618	1	
50 cm		5 mm	V1103620	1	
60 cm		6 mm	V1103626	1	
		7 mm	V1103627	1	
		8 mm	V1103628	1	
70 cm		5 mm	V1103630	1	
80 cm		6 mm	V1103636	1	
		7 mm	V1103637	1	
		8 mm	V1103638	1	
Thin wall, helix		50 cm	5 mm	V1103670	1
		60 cm	6 mm	V1103676	1
	7 mm		V1103677	1	
	8 mm		V1103678	1	
	70 cm	5 mm	V1103680	1	
	80 cm	6 mm	V1103686	1	
		7 mm	V1103687	1	
8 mm		V1103688	1		